

Core body temperature monitoring during daily life

greenTEG AG

2020



Study

Core body temperature monitoring during daily life

Goal: Detection daily temperature cycles using greenTEG's smartwatch integrated sensor

Activity: Free-living

Requirements: 2x 1h of high-intensity sports


Position: Wrist

Reference: Ingestible radio pill

Data collected: > 2.5 million data points

Calibration: No calibration required

First reading: After 4 minutes





Working



Eating







3

Sleeping

Sports



Results:

Algorithm 1: 1 XU integrated into Health Watch

Sensors used:

- 1 XU greenTEG's Thermal Energy Transfer Sensor
- 3rd Party PPG sensor

Results:


- High accuracy of core body temperature monitoring on the wrist
- Capable of tracking the circadian cycle
- Good compensation of any thermal influences in the environment

Statistics (over all measurements):

- Mean absolute deviation: 0.30 °C
- Standard deviation: 0.26 °C
- Correlation Coefficient: 0.75

Calibration: No calibration required

First reading: After 4 minutes





Results:

Algorithm 2: greenTEG's Demonstration Watch

Sensors used:

2 XU - greenTEG's Thermal Energy Transfer Sensors

Results:

- High accuracy of core body temperature monitoring on the wrist
- Capable of tracking the circadian cycle
- Good compensation of any thermal influences in the environment

Statistics (over all measurements):

Mean absolute deviation: 0.22 °C

Standard deviation: 0.29 °C

Correlation Coefficient: 0.76

Calibration: No calibration required

First reading: After 4 minutes

